

# **MATERIAL SAFETY DATA SHEET**

## **Frozen Cell Cultures**

#### Material Safety Data Sheet for:

UK Health Security Agency (UKHSA) frozen cell cultures

Issued to: Users of UKHSA frozen cell cultures

Access: Document to be downloaded from Culture Collections website https://www.culturecollections.org.uk

Culture Collections UK Health Security Agency Porton Down Salisbury SP4 0JG UK Telephone: +44 (0) 1980 612512 https://www.culturecollections.org.uk/contactus

Page 1 of 7



### MATERIAL SAFETY DATA SHEET FOR UKHSA FROZEN CELL CULTURES

#### 1. Identification of the substance/mixture and of the company/undertaking

**Product name:** Various human and animal cell cultures supplied by the European Collection of Authenticated Cell Cultures (ECACC) which is part of the UK Health Security Agency Culture Collections. This includes a collection of Induced pluripotent stem cells (iPSCs) which have been produced using genetic modification. Please consult your local health and safety guidelines to ensure you can safely handle these products.

Refer to the relevant cell line data entry on the Culture Collections website: <u>https://www.culturecollections.org.uk</u>

Contact: Culture Collections UK Health Security Agency Porton Down Salisbury SP4 0JG, UK Telephone (working hours): +44 (0)1980 612512 Telephone (out of hours): +44 (0)1980 612100 https://www.culturecollections.org.uk/contactus

#### 2. Hazards Identification

#### **Chemical Hazards:**

Frozen cultures may contain 5 to 10% (v/v) dimethyl sulphoxide (DMSO). DMSO may be harmful and toxic if in contact with skin or ingested, (R23/24/25). It also maybe irritating to eyes and respiratory system (R36/37/38). Thawed contents of vials should not come into contact with skin, eyes or digestive and respiratory epithelium (S24/25) and should be diluted upon use with culture media. Persons handling vials of frozen cells containing DMSO should wear a laboratory overall, protective glasses and insulated gloves (S36/37).

#### **Biological hazards:**

Human or animal cell lines may not be known to contain specific human pathogens but the possibility of the presence of a pathogenic contaminant or pathogenic adventitious virus in the cells cannot be excluded. Therefore, it is recommended that all human and animal cell lines are handled at containment Level 2 as Advisory Committee on Dangerous Pathogens (ACDP) Hazard Group 2 organisms, unless a higher ACDP Hazard Group is specified.

The relevant cell line data entry on the Culture Collections website includes any specific instructions applicable to the biohazardous potential of this cell line. These

Page 2 of 7



instructions should be considered by the user when assessing risks of the use of the cells. Users should refer to the references in the relevant cell line data entry on the Culture Collections website: <u>https://www.culturecollections.org.uk</u>. These cell lines have not been screened for adventitious agents.

#### Health Effects:

- Eyes: Unknown
- Skin: Unknown
- Ingestion: Unknown
- Inhalation: Unknown

#### **Physical Hazards:**

Where cell lines are shipped as frozen vials there is a small risk that the vial may be pressurised, due to the expansion of trapped liquid nitrogen and could explode on warming. Such a risk will be increased if the vial has been shipped to the customer in a liquid nitrogen container (dry-shipper). It is recommended that persons handling vials of frozen cells should wear a laboratory overall, protective glasses and protective laboratory gloves.

This sheet does not constitute an assessment as required by the Control of Substances Hazardous to Health Regulations 1994.

The information contained in this publication is given in good faith and is accurate to the best of our knowledge.

#### 3. Composition/information on ingredients

Frozen cultures contain 5 to 10% (v/v) dimethyl sulphoxide (DMSO) the remainder consisting of Foetal Bovine Serum (FBS). Human or animal cell lines may not be known to contain specific human pathogens but the possibility of the presence of a pathogenic contaminant or pathogenic adventitious virus in the cells cannot be excluded. Therefore, it is recommended that all human and animal cell lines are handled at containment Level 2. as ACDP Hazard Group 2 organisms unless a higher ACDP Hazard Group is specified.

#### 4. First aid measures

If accidental contact with material occurs laboratory staff must follow the local first aid procedures that are normally applied following exposure to organisms of ACDP Hazard Group 2.

- Eyes: Irrigate with physiological saline or water and seek medical advice immediately
- Skin: Wash thoroughly with soap and water and seek medical advice immediately
- Ingestion: Seek medical advice immediately
- Inhalation: Seek medical advice immediately

Page **3** of **7** 



#### 5. Fire-fighting measures

- Extinguisher medium: Use medium suitable for surrounding environment
- Unsuitable extinguisher medium: N/A
- Protective equipment for firefighting: N/A

#### 6. Accidental release measures

Personal precautions: avoid direct contact with the thawed material. Do not open the primary containers unless authorised to do so. Wear a laboratory overall, protective laboratory gloves and safety glasses.

Environmental precautions: if spillage occurs place absorbent material over the spillage and disinfect. See below.

Spillage of thawed material: wear a laboratory coat, safety glasses and protective laboratory gloves. Place paper towels or other absorbent material over the spillage. Pour disinfectant over spillage to saturate and leave for 30 minutes prior to cleaning and disposal. The preferred disinfectant is 10% v/v sodium hypochlorite (10,000 parts per million available chlorine). This should not be used in combination with other disinfectants, see your local risk assessment or contact the manufacturer of the disinfectant for additional information.

#### 7. Handling and storage

Personal protective equipment comprised of laboratory coat, protective laboratory gloves and safety glasses should be worn when handling (unpacking) human and animal cell lines. The dry ice (solid carbon dioxide) used to ship frozen vials should be allowed to evaporate in a well-ventilated area. Do not dispose of dry ice in a sealed container.

Vials or flasks containing human and animal cells should be opened in a Class II microbiological safety cabinet under conditions of ACDP Hazard Group 2.

Induced pluripotent stem cells have been produced using genetic modification please consult your local health and safety guidelines to ensure you are able to safely handle this product.

Detailed discussions of laboratory safety procedures are provided in: "Laboratory Safety: Principles and Practice" (Fleming, et al, 1995); the Journal of Tissue Culture Methods (Caputo, 1988), and in the U.S. Government Publication, "Biosafety in Microbiological and Biomedical Laboratories (BMBL) 5<sup>th</sup> Edition" (CDC, 2009). This

Page 4 of 7



publication is available on the Center for Disease Control, Office of Health and Safety's website: www.cdc.gov/biosafety/publications/bmbl5/

#### 8. Exposure controls / personal protection

Engineering control measures: Vials containing human and animal cells should be opened in a Class II microbiological safety cabinet under conditions of ACDP Hazard Group 2. Personal protective equipment comprised of laboratory coat, protective laboratory gloves and safety glasses should be worn.

Respiratory protection: avoid aerosol production and inhalation. Handle as for ACDP Hazard Group 2.

Hand protection: always wear protective laboratory gloves

Eye protection: always wear safety glasses

#### 9. Physical and chemical properties

- Appearance: Frozen fluid in small glass or plastic containers (vials). The majority of cell cultures are supplied in plastic vials
- Solid/liquid/gas: Solid (frozen state)
- Appearance: yellow or pink solid for frozen cultures. Aqueous pH 6-8
- The product is provided as a frozen culture of animal cells
- The frozen components may include but are not limited to water, inorganic salts, vitamins, amino acids, carbohydrates, lipids, proteins (animal-derived) and cryoprotectant (dimethyl sulphoxide 10% v/v; or glycerol 10-20% v/v), phenol red.

#### 10. Stability and reactivity

Reactivity data: Stable. Hazardous polymerization will not occur.

#### 11. Toxicological information

- Routes of exposure: Not applicable
- Acute effects: Not applicable
- Chronic effects: Not applicable
- Special considerations: In its thawed liquid state this material is not normally toxic but avoid aerosol formation and inhalation. Vials contain dimethyl sulphoxide 10% v/v which is an irritant that readily penetrates the skin.

#### 12. Ecological information

Page 5 of 7



- Mobility: May apply in certain cases of genetic modification consult the relevant cell line data entry on Culture Collections website: <u>https://www.culturecollections.org.uk</u>
- Persistence / degradability: N/A
- Bioaccumulation: N/A
- Ecotoxicity: N/A

#### 13. Disposal considerations

Follow established procedures for Containment (Biosafety) Level 2 and follow all national, regional and local regulations. The UK Environmental Protection Act 1990 applies.

#### Methods for disposal for thawed content:

- Spillage: wear a laboratory coat, safety glasses and protective laboratory gloves. Place paper towels or other absorbent material over the spillage. Pour disinfectant over spillage to saturate and leave for 30 minutes prior to cleaning and disposal. The most appropriate disinfectant is 10% v/v sodium hypochlorite (10,000 parts per million available chlorine). This should not be used in combination with other disinfectants. Check your local risk assessment or contact the manufacturer of the disinfectant for additional information
- Waste disposal: Decontaminate prior to disposal with a 10% sodium hypochlorite solution and dispose of decontaminated liquid waste down a designated sink with running water. Solid waste should be placed in a sealed bag and labelled and destroyed by incineration.

#### 14. Transport information

Additional information arising from the Carriage of Dangerous Goods by Road & Air (Classification, Packaging and Labelling) Regulations:

UN 1845: Dry Ice. Dry ice is not subject to ADR (road regulations) except for 5.5.3. It is only classified as a dangerous goods under IATA regulations (air)

Where a dry shipper is used containing liquid nitrogen fully absorbed in a porous material IATA (International Air Transport Association) Dangerous Goods Regulations do not apply.

Most ECACC cell lines are not classified as dangerous goods as they are considered non-infectious to humans or animals and are not genetically modified; therefore, they are not subject to IATA or European Agreement concerning the International Carriage of Dangerous Goods by Road (ADR) regulation for dangerous goods.

When the following categories apply, ECACC will ensure the outer packaging indicates the appropriate packaging requirements:

#### Page 6 of 7



- Biological Substance Category B UN3373 packed in compliance with IATA packing instruction 650
- Genetically Modified Organisms (GMOs) UN3245 packed in compliance with IATA packing instruction 959

#### 15. <u>Regulatory information</u>

ECACC confirms that all necessary licences (import, holding, transfer and export) required for the consignment of this material are in place. The recipient is only required to provide evidence of permits and licences to receive and handle ACDP Hazard Group 3 restricted materials.

For the relevant licences consult the Culture Collections website: <u>https://www.culturecollections.org.uk</u>.

This organism/material may be covered by United Kingdom or International legislation.

#### 16. Other information

In the event of an accident involving exposure of a person to the material contained in the samples, contact the Culture Collections (+44 (0)1980 612512) during normal UK working hours. Refer to section 1 for full contact details.

The above information is correct to the best of our knowledge. All materials and mixtures may present unknown hazards and should be used with caution.

The user should make independent assessments and decisions regarding the completeness of the information based on all sources available.

The Culture Collections shall not be held liable for any damage resulting from handling or contact with the above product.

Page 7 of 7